

<b>Title</b>	LP20E Service BOM		<b>Page</b>	1 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

This document is electronically signed in the Physio-Control, Inc. Product Lifecycle Management (PLM) System. Approvals can be obtained from the PLM system which displays the selected approvers, their approval roles, and approval dates.

**PROPRIETARY AND CONFIDENTIAL DOCUMENT FOR PHYSIO-CONTROL, INC. USE ONLY.** This document is property of Physio-Control and may not be used, reproduced, published or disclosed to others without authorization from Physio-Control.



Physio-Control, Inc.  
 11811 Willows Road NE, P.O. Box 97006, Redmond, WA 98052  
 Tel 425.867.4000 Toll-free 800.442.1142  
[www.physio-control.com](http://www.physio-control.com)

Owner

Corporate

<b>Attachments (Optional):</b>	<b>Qty</b>	<b>Pages</b>	<b>Page Sequencing</b>
3343496_AC_LP20E_Service_BOM.pdf	1	18	1-18

**Change Description**

BOM Changes:

- Updated dash number for Service Manual part number (3314176-009) in section 1 and 2.
- Reference designators updated to align with current designators used in Service Manual.
- Numerical reference designators removed and replaced with alphanumeric designators to align with Service Manual.
- Updated several outdated part descriptions to match description of current revision of part.
- Updated catalog number and dash number for 3337495-003.
- Removed A11 (3317230-002) as it is no longer used.

New part numbers added to service BOM:

3201374-012, 3202018-000, 3201415-004, 3202718-035, 3202718-036, 3202718-037, 3202718-038, 3202718-039, 3202718-040, 3202718-044, 3202718-045, 3202718-046, 3202718-047, 3202718-048, 3202718-054, 3202718-055, 3202718-056, 3308907-000, 202253-761, 3317227-003, 3317230-002, 3201643-033, 3201643-026, 3205278-001, 3309555-000, 3206469-003, 3200996-002

Parts were updated in 3343496-000 Rev AC EBOM to align with changes made in this revision of the LIFEPAK 20e Service BOM. Refer to redlines for details of change. 3309555-000, 3308907-000, 202253-761, and 3201374-012 were previously connected to the e-BOM at Rev A but not present in this service BOM document at Rev A. All other part numbers were added to both the e-BOM and service BOM document at Rev AC.

- Further revisions were needed to drawing 3343496 Rev AB. Rev AB will be obsoleted on CO 1312439 and Rev AC will be released. 3343496-000 Rev AC was also created and attached to 3343496 Rev AC.

Note: This document along with other service-related documents are being switched to double-letter revisions due to a need associated with the training software that the service department uses.

<b>Title</b>	LP20E Service BOM		<b>Page</b>	2 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

1. **PURPOSE**

This document provides the list of service parts for the LIFEPAK 20E monitor/defibrillator and CodeManagement Module and provides cross references between reference designators (ex. A01) used in the LIFEPAK 20E Service Manual (3314176-009), and part and catalog numbers.

2. **SCOPE**

Scope is limited to individual parts and service repair kits which are available for service purposes for the LP20E Service Manual (3314176-009). Kit contents are not listed in this document and exist in electronic bills of materials (eBOMs) in the Product Lifecycle Management (PLM) system.

3. **SERVICE BOM PART**

Item Number	Description
3343496-000	Bill of Material, Service, LIFEPAK 20E

See 3343496-000 EBOM for parts list or reference the tables in Sections 4 and 5.

4. **SERVICE INDIVIDUAL PARTS LIST**

Note: Reference Designators are unique identifiers for the part used to cross-reference between the Service Manual or Service Process Work Instruction and the Service BOM. The reference designators have been assigned based on the following prefix system:

- A prefix- Printed Circuit Board Assemblies (PCBAs), keypads, and modules
- W prefix – Cables, subassemblies
- F prefix – Fasteners (e.g. screws, nuts, and standoffs)
- M prefix – Plastic parts and mechanical components
- E prefix – Electrical components
- L prefix – Labels and adhesives
- K prefix – Kits, or parts assembled into a kit for a specific repair.

Note: For some reference designators, there are multiple language options (e.g. labels and keypads). Refer to the tables below to select the version that applies.

<b>Title</b>	LP20E Service BOM			<b>Page</b>	3 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b>	AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023		

Click on a link to navigate to a parts list:

[-Top Case Parts List](#)

[-AED Door/Latch Label Kits](#)

[-Front Case Parts List](#)

[-A15 Elastomer Keypad- Languages](#)

[-Boardstack Parts List](#)

[-Standard Paddles Parts List](#)

[-Bottom Case Parts List](#)

[-Standard Paddles Label Languages](#)

[-Device Labeling Including Label Set \(12\)](#)

[-Charge Button Languages](#)

[-LP20E Device Label Set Languages](#)

[-CMM Replacement Parts List](#)

[-Manual Latch Label Languages](#)

[-LP20E Service Repair Kits](#)

[-Adult VF Dose Label Languages](#)

[-AED Door/Latch Label Kits](#)

Top Case Parts List				
Reference Designator	Part Number	Catalog Number	Description	Notes
F01	202253-592	21300-005334	SCREW,MACHINE,PANHEAD,NY LOK,6-32 X 1.750	
F02	3206185-001	21300-001032	SCREW-PH,SS,P,NYLOK,6-32 X .375	
M01	3202497-002	21330-001036	CASE-TOP,SUB ASSEMBLY,LP20	

Front Case Parts List				
Reference Designator	Part Number	Catalog Number	Description	Notes
A06	3337495-003	21330-001642	PCB ASSY - PROGRAMMED USER INTERFACE, LP20	Part of kit REF 40402-000037, 40402-000038, 40402-000039
E01	3202033-500	21300-006680	INV-BKLIGHT,CCFT,LCD,1.5KV OUT,8TO20VIN,ROHS	Part used with E03 - REF 21300-007363 Part of kit REF 40402-000010
A08	3308907-000	21300-008166	MODULE-OEM,LED DRIVER	Part used with E04 - REF 21300-008077 Part of kit REF 40402-000027
E03	3205278-001	21300-007363	DISPLAY, COLOR LCD, KYOC ERA, 5.7 , ROHS	Part of kit REF 40402-000009

<b>Title</b>	LP20E Service BOM		<b>Page</b>	4 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

Front Case Parts List				
E04	3309555-000	21300-008077	LCD DISPLAY - COLOR , 5. 7IN , ROHS	Requires REF 21330-001345 (RoHS) or greater and Redux Device S/W Part of kit REF 40402-000027
A07	Multiple	Multiple	Elastomer Keypad	Refer to tables for multiple keypad language options.
E05	3319086-001	21300-008133	SPEAKER ASSEMBLY,ROHS	RoHS
W04	3201145-000	21300-004264	ENCODER-OPTICAL, ROTARY, 16 POSN, 3.3V, 4IN, ROHS	
W15	3300152-000	21300-007774	CABLE-ASSEMBLY, ACTIVE COLOR, ROHS	
W16	3200996-002	21300-004237	CABLE-BACKLIGHT INVERTER ,ACTIVE COLOR	
W17	3201000-004	21300-008056	FLEX ASSEMBLY - UI TO STACK	
M02	3200624-009	21300-007086	CASE-FRONT,PAD PRINTED,LP20	Part of kit REF 40402-000002
L01	Multiple	Multiple	Label set (6 labels)	Refer to tables for multiple language options.
M03	3200633-000	21300-004620	KNOB-SPEED DIAL,LP20	
M04	3200638-001	21300-004837	COVER PLATE-DOOR BLANK ASSY WITH ADHESIVE TAPE	Part of kit REF 40402-000002
L02	3201499-900	21501-000767	LABEL-MANUAL LATCH,ENG,SPA,POR,BRA	Select other language
M05	3200637-001	21300-004836	PLATE-DOOR HINGE ASSY,WITH ADHESIVE TAPE	Part of kit REF 40402-000002
M06	3201610-000	21300-004649	BRACKET-MOUNTING,SPEAKER	
M07	3200913-003	21300-004233	LENS - DISPLAY, LIFEPAK 20	Part of various kits
M08	3200640-015	21300-004838	BRACKET - DISPLAY, LP20, ACTIVE, ADHESIVE / FOAM	Part of various kits
F03	200805-000	21300-000584	NUT-HEX,SS,LOCK 4-40 X.250W	
F04	202253-761	21300-001038	SCREW-M,CS,Z,PH, NYLOCK, 4-40 X .312L	
M09	3202056-004	21300-004252	DOOR ASSEMBLY - AED, LIFEPAK 20	Pat of kit MIN 3202360. See AED Door/Latch Label Kits
L13	3201111-003	21300-004241	PAD-ELASTOMER,THERMALLY CONDUCTIVE,BACKLIGHT INV	
M10	3202246-003	21300-004254	HARNESS-GROUNDING STRAP	Active Display
M11	3202246-004	21300-004255	HARNESS-GROUNDING STRAP,	User Interface PCB

<b>Title</b>	LP20E Service BOM		<b>Page</b>	5 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

<b>Front Case Parts List</b>				
			MODIFIED	
L14	804447-041	21300-004807	SPACER-FOAM,LP20 UI PCB	Part of various kits
M12	3202246-005	21300-004884	HARNESS-GROUNDING STRAP, SPEED DIAL	
F05	3205497-308	21300-006141	RVT-NYL,SNAP,0.118 DIA,0 .31LG,AEH 06-01-3080	Part of various display kits
E15	3202018-000	21300-004250	CABLE-EXTENDER,ACTIVE DISPLAY	
L03	Multiple	Multiple	Label - Adult VF Dose	Refer to tables for multiple language options

<b>Title</b>	LP20E Service BOM		<b>Page</b>	6 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

<b>Boardstack Parts List</b>				
<b>Reference Designator</b>	<b>Part Number</b>	<b>Catalog Number</b>	<b>Description</b>	<b>Notes</b>
A01	3337492-002	21330-001561	PCB ASSY-PROGRAMMED SYSTEM CONTROLLER, LP20E	Part of kit REF 40402-000037 and REF 40402-000040
A02	3337384-001	21330-001558	PCB ASSY - PP, 3 CHANNEL WITH SIR, LP20-LP20E	Part of kit REF 40402-000045 (4th Edition)
A03	3337489-100	21330-001560	PCB ASSY - THERAPY MODULE, LP20	Part of kit REF 40402-000028 (RoHS)
A03	3337489-200	21330-001563	PCB ASSY - THERAPY, SIM MOD NOT INSTALLED, LP20	Part of kit REF 40402-000029 (RoHS)
A04	3337071-001	21330-001557	PCB ASSY-OEM INTERFACE, MASIMO	Part of kit REF 40402-000044 (4th Edition)
A05	3318279-001	21996-000074	MODULE-OEM, PULSE OXIMETER, 3780, ROHS	Part of kit REF 40402-000034 (RoHS)
A05	3318279-002	21996-000117	MODULE-OEM, PULSE OXIMETER WITH MNC, 3783, ROHS	For devices with Masimo Nellcor compatibility, Order kit REF 40402-000035 (RoHS)
E06	3202773-000	21300-005068	INDUCTIVE RESISTOR-ASSEMBLY,LP20	
W02	3202383-001	21300-004307	WIRE-CAPACITOR DISCHARGE	
M13	3201374-010	21300-004242	MOUNT-ISOLATION,PCB, THERAPY	
F06	202253-772	21300-006430	SCREW,MACHINE,PANHEAD,NY LOK,4-40 X 1.000	
M14	3200927-006	21300-004236	SHIELD-PCB STACK	Part of kit REF 40402-000017
M15	3201375-006	21300-004815	STANDOFF - MALE/MALE, .250 HEX, .375L	Part of kit REF 40402-000017
W07	3201007-002	21300-004704	CONNECTOR-BOARD STACKER, MODIFIED	Part of kit REF 40402-000033 (RoHS)
M16	3200639-007	21300-004228	BRACKET,SUPPORT,PCB	
M17	3201375-005	21300-004816	STANDOFF - MALE/MALE, .250 HEX, .250L	
L15	3201111-000	21300-007457	PAD-ELASTOMER,THERMALLY CONDUCTIVE,THERAPY	No longer used in assembly
M18	3201374-012	21300-004243	MOUNT-ISOLATION,PCB,PARAMETER	Part of kit REF 40402-000017
F03	200805-000	21300-000584	NUT-HEX,SS,LOCK 4-40 X.250W	

<b>Title</b>	LP20E Service BOM		<b>Page</b>	7 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

<b>Boardstack Parts List</b>				
F04	202253-761	21300-001038	SCREW-M,CS,Z,PH, NYLOCK, 4-40 X .312L	
F07	202253-764	21300-004599	SCREW,MACHINE,PANHEAD,NYLOK,4-40 X .500	
M19	200536-001	21300-000499	RETAINR-CABLE TIE,NYLON, .10 W, X 4.0 L	
F08	201874-002	21300-008019	SCREW-M,PH,CP,P,SEMS,#4- 40,5-16 L,SPLIT WASHER	
F09	3202489-031	21300-005120	SCREW-CAP, SKT HD, HEX REC, 4-40, 0.312L, NYL	3/32 hex drive
M20	3203897-000	21300-006038	SHIELD - NOMEX	Part of kit REF 40402-000017
F10	200804-102	21300-000580	WASHER,FLAT,ROUND,S,ZINC ,.312D,.125ID	
M21	3201374-009	21300-005578	MOUNT-ISOLATION,PCB,SYSTEM CONTROLLER	Part of kit REF 40402-000033 (RoHS)
M22	3206405-000	21300-006593	SHIELD - EMI, PCB STACK	Part of kit REF 40402-000017
M42	3201374-012	21300-005187	MOUNT-ISLOATION,PCB,PARAMETER	
M43	3201415-004	21300-004245	BRACKET-INDUCTIVE RESISTOR,LP20	
A10	3317227-003	21330-001354	PCB ASSY-PROGRAMMED SYSTEM CONTROLLER,LP20	
E07	202305-000	21300-001052	BATTERY-LITH,3V,COIN	

<b>Title</b>	LP20E Service BOM		<b>Page</b>	8 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

Bottom Case Parts List				
Reference Designator	Part Number	Catalog Number	Description	Notes
E08	3317740-002	21330-001569	POWER SUPPLY ASSEMBLY-LIFEPAK 20E,4TH ED,NEW MCU	Part of Kit REF 40402-000052
E09	3205296-006	11141-000112	BATTERY PACK - LI-ION, LP20E, E-CELL	
A12	3318747-001	21240-000001	PRINTER-50MM, THERMAL PRINTHEAD, ROHS	
E10	3200846-002	21300-004232	CAPACITOR-ENERGY STORAGE , 2.25 IN DIA	Part of kit REF 40402-000021
W01	3200474-009	21300-007366	CONNECTOR ASSY-THERAPY RECPT,HVY MTL LP20E,ROHS	Part of kit REF 40402-000036 (RoHS)
W03	3200926-007	21300-004235	FLEX ASSEMBLY - IRDA,ROHS	(RoHS)
W05	3200925-006	21300-004234	FLEX ASSY - SPO2	Part of Masimo SpO2 Connector Repair Kit REF 40402-000014
W06	3318315-000	21330-001521	RECEPTACLE - ECG, 7 CONTACT, WIRED, ROHS	(RoHS)
W08	3206579-001	21330-001166	CABLE ASSEMBLY-BATTERY INTERFACE	Part of kit REF 40402-000032 (RoHS)
W09	3201241-000	21300-004669	CABLE-POWER PCBA-26 PIN TO THERAPY PCBA	
W10	3206857-000	21330-001165	CABLE ASSY-POWER/THERAPY	Part of kit REF 40402-000032 (RoHS)
W11	3201997-506	21300-007072	CABLE ASSY- D-SUB CONN, LP20/LP20E,ROHS	Part of kit REF 40402-000032 (RoHS)
W12	3206469-003	21330-001164	CABLE ASSY-OUTPUT,POWER SUPPLY	Part of kit REF 40402-000032 (RoHS)
W13	3201001-006	21300-008058	FLEX ASSEMBLY - PRINTER	
W14	3319086-000	21300-008132	SPEAKER ASSEMBLY EXTENSION,ROHS	(RoHS)
L01	Multiple	Multiple	Label Set (6 labels), LIFEPAK 20e	Refer to tables for multiple language options
F11	3202228-562	21300-003883	STDF-HEX,THD,M/F,4-40,0.188DIA,562L,CS,ZN,NYLOCK	Part of kit REF 40402-000020
M23	3201408-000	21300-004621	BEZEL-PRINTER,LP20	
M24	3200626-004	21300-004306	SHROUD-PRINTER,LP20	
L16	804447-36	21300-007458	SPACER - FOAM,CAPACITOR BKTS	
M25	3200922-000	21300-004619	BRACKET-CAPACITOR SUPPORT	Part of kit REF 40402-000020
M26	3201597-000	21300-004653	GASKET-EMI,FOAM CORE	
M27	3200628-004	21300-006412	DOOR-LI-ION BATTERY WITH FOOT	



<b>Title</b>	LP20E Service BOM		<b>Page</b>	9 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

<b>Bottom Case Parts List</b>				
			ASSY, LP20E	
M28	802885-00	21300-002137	FOOT-MOUNTING,LP8/9/10	Part of bottom case assembly
M29	200040-001	21300-000149	RING-RETAINING,CRES,1.125IN SHAFT DIA	Part of kit REF 40402-000036 (RoHS)
M30	3200921-000	21300-004602	CLIP-CONNECTOR MOUNTING,SPO2	Part of kit REF 40402-000014
M31	3200625-006	21300-004889	CASE-BOTTOM,FOOTED,LP20	
L17	3009787-003	21300-004110	SHIELD - DIELECTRIC, CAPACITOR	Part of REF 40402-000021. Part of E10 Assembly.
F04	202253-761	21300-001038	SCREW-M,CS,Z,PH, NYLOCK, 4-40 X .312L	
L04	3201643-025	21501-002888	LABEL-SERIAL NUMBER, LP20E	Refer to Labels Assembly
M11	3202246-004	21300-004255	HARNES-GROUNDING STRAP, MODIFIED	
F12	3202377-010	21300-004400	FSTNR-DRIVE,SNAP-IN,0.126DIA HOLE,ACETAL	
M41	802885-01	21300-002138	FOOT-MOUNTING	Attached to battery door (103)
M19	200536-001	21300-000499	RETAINR-CABLE TIE,NYLON, .10 W, X 4.0 L	Part of kit REF 40402-000021
M32	3206900-001	21300-006962	SHIELD-NOMEX,DEFIB CAPACITOR	Part of kit REF 40402-000020
M33	3300694-000	21300-007794	PORON FOAM - LP20,LP20E CAPACITOR	Part of kit REF 40402-000021
M34	3203445-001	21300-005784	SEAL, THERAPY CONNECTOR MOUNT, LP20	Part of kit REF 40402-000036 (RoHS)
L05	3206152-000	21501-001625	LABEL-NIMH BATTERY WARNING, LP20E	

<b>Title</b>	LP20E Service BOM		<b>Page</b>	10 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

Device Labeling Including Label Set (12)				
Reference Designator	Part Number	Catalog Number	Description	Notes
L06	3201290-004	21501-001029	LABEL-REAR CONNECTOR,LP20/LP20E	A in Labels Assembly Diagram
L07	3201274-030	21501-000894	LABEL-FRONT CONNECTOR, MASIMO, LP20, ENGLISH	B in Labels Assembly Diagram
L07	3201274-015	21501-000879	LABEL-FRONT CONNECTOR, NO SPO2, LP20, ENGLISH	B in Labels Assembly Diagram
L08	3201275-003	21501-000909	LABEL-PHYSIO-CONTROL LOGO LP20 LP20E	C in Labels Assembly Diagram
L08	3201273-016	21501-001622	LABEL-LOGO,LP20E, ENGLISH	D in Labels Assembly Diagram
L09	3009060-001	21501-000088	LABEL-FDA	E in Labels Assembly Diagram
L10	3202375-003	21501-000779	LABEL-MASIMO PATENT	F in Labels Assembly Diagram
L02	Multiple	Multiple	Manual Latch Label	Refer to tables for multiple language options
L04	3201643-025	21501-002888	LABEL-SERIAL NUMBER, LP20E	SN label illustration
L04	3201643-033	21501-002992	LABEL-SERIAL NUMBER, LP20E, 3RD ED, NEW FORMAT	SN label illustration
L04	3201643-026	21501-002889	LABEL-SERIAL NUMBER, LP20E, NO CE	SN label illustration
L03	3206926-001	21501-001935	LABEL-ADULT VF DOSE, AHA , ENGLISH	Select language
L11	3320539-000	21501-002882	LABEL - UDI	See Labels Assembly diagram

LIFEPAK 20E Device Label Set Languages				
Reference Designator	Part Number	Catalog Number	Description	Notes
L01	3206034-008	21501-002660	LABEL SET - LP20E, ENG, 3RD EDITION	
L01	3206034-009	21501-002661	LABEL SET - SPO2, LP20E, ENG, 3RD EDITION	
L01	3206034-086	21501-002666	LABEL SET - LP20E, FRE, 3RD EDITION	
L01	3206034-087	21501-002667	LABEL SET - SPO2, LP20E, FRE, 3RD EDITION	
L01	3206034-046	21501-002662	LABEL SET - LP20E, GER, 3RD EDITION	
L01	3206034-047	21501-002663	LABEL SET - SPO2, LP20E, GER, 3RD EDITION	
L01	3206034-126	21501-002670	LABEL SET - LP20E, SPA, 3RD EDITION	
L01	3206034-127	21501-002671	LABEL SET - SPO2, LP20E, SPA, 3RD EDITION	
L01	3206034-066	21501-002664	LABEL SET - LP20E, ITA, 3RD EDITION	
L01	3206034-067	21501-002665	LABEL SET - SPO2, LP20E, ITA, 3RD EDITION	

<b>Title</b>	LP20E Service BOM			<b>Page</b>	11 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b>	AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023		

LIFEPAK 20E Device Label Set Languages				
L01	3206034-186	21501-002676	LABEL SET - LP20E, SWE, 3RD EDITION	
L01	3206034-187	21501-002677	LABEL SET - SPO2, LP20E, SWE, 3RD EDITION	
L01	3206034-206	21501-002678	LABEL SET - LP20E, DAN, 3RD EDITION	
L01	3206034-207	21501-002679	LABEL SET - SPO2, LP20E, DAN, 3RD EDITION	
L01	3206034-106	21501-002668	LABEL SET - LP20E, DUT, 3RD EDITION	
L01	3206034-107	21501-002669	LABEL SET - SPO2, LP20E, DUT, 3RD EDITION	
L01	3206034-226	21501-002680	LABEL SET - LP20E, FIN, 3RD EDITION	
L01	3206034-227	21501-002681	LABEL SET - SPO2, LP20E, FIN, 3RD EDITION	
L01	3206034-246	21501-002682	LABEL SET - LP20E, NOR, 3RD EDITION	
L01	3206034-247	21501-002683	LABEL SET - SPO2, LP20E, NOR, 3RD EDITION	
L01	3206034-266	21501-002684	LABEL SET - LP20E, POL, 3RD EDITION	
L01	3206034-267	21501-002685	LABEL SET - SPO2, LP20E, POL, 3RD EDITION	
L01	3206034-146	21501-002672	LABEL SET - LP20E, POR, 3RD EDITION	
L01	3206034-147	21501-002673	LABEL SET - SPO2, LP20E, POR, 3RD EDITION	
L01	3206034-166	21501-002674	LABEL SET - LP20E, BRA, 3RD EDITION	
L01	3206034-167	21501-002675	LABEL SET - SPO2, LP20E, BRA, 3RD EDITION	
L01	3206034-406	21501-002696	LABEL SET - LP20E, JPN, 3RD EDITION	
L01	3206034-407	21501-002697	LABEL SET - SPO2, LP20E, JPN, 3RD EDITION	
L01	3206034-346	21501-002692	LABEL SET-LP20E, CHINESE, 3RD EDITION	
L01	3206034-347	21501-002693	LABEL SET-SPO2, LP20E, CHINESE, 3RD EDITION	
L01	3206034-286	21501-002686	LABEL SET - LP20E, HUN, 3RD EDITION	
L01	3206034-287	21501-002687	LABEL SET - SPO2, LP20E, HUN, 3RD EDITION	
L01	3206034-306	21501-002688	LABEL SET - LP20E, CZE, 3RD EDITION	
L01	3206034-307	21501-002689	LABEL SET - SPO2, LP20E, CZE, 3RD EDITION	
L01	3206034-326	21501-002690	LABEL SET - LP20E, RUS, 3RD EDITION	
L01	3206034-327	21501-002691	LABEL SET - SPO2, LP20E, RUS, 3RD EDITION	
L01	3206034-386	21501-002694	LABEL SET-LP20E, KOREAN, 3RD EDITION	
L01	3206034-387	21501-002695	LABEL SET-SPO2,LP20E, KOREAN, 3RD EDITION	

Manual Latch Label Languages				
Reference Designator	Part Number	Catalog Number	Description	Notes
L02	3201499-900	21501-000767	LABEL-MANUAL LATCH,ENG,SPA,POR,BRA	
L02	3201499-902	21501-000924	LABEL-MANUAL LATCH,FRE,DAN	
L02	3201499-901	21501-000925	LABEL-MANUAL LATCH,GER,SWE,NOR	
L02	3201499-060	21501-000932	LABEL-MANUAL LATCH,ITALIAN	
L02	3201499-100	21501-000930	LABEL-MANUAL LATCH,DUTCH	
L02	3201499-220	21501-000929	LABEL-MANUAL LATCH,FINNISH	
L02	3201499-260	21501-000928	LABEL-MANUAL LATCH,POLISH	

<b>Title</b>	LP20E Service BOM		<b>Page</b>	12 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

Manual Latch Label Languages				
Reference Designator	Part Number	Catalog Number	Description	Notes
L02	3201499-400	21501-000926	LABEL-MANUAL LATCH,JAPANESE	
L02	3201499-340	21501-000927	LABEL-MANUAL LATCH,CHINESE	
L02	3201499-280	21501-001349	LABEL-MANUAL LATCH, HUNGARIAN	
L02	3201499-300	21501-001350	LABEL-MANUAL LATCH, CZECH	
L02	3201499-320	21501-001351	LABEL-MANUAL LATCH, RUSSIAN	
L02	3201499-380	21501-001352	LABEL-MANUAL LATCH, KOREAN	

Adult VF Dose Label Languages				
Reference Designator	Part Number	Catalog Number	Description	Notes
L03	3206926-001	21501-001935	LABEL-ADULT VF DOSE, AHA , ENGLISH	
L03	3206926-081	21501-001858	LABEL-ADULT VF DOSE, AHA 05, FRENCH	
L03	3206926-041	21501-001856	LABEL-ADULT VF DOSE, AHA 2005,GERMAN	
L03	3206926-121	21501-001860	LABEL-ADULT VF DOSE, AHA 2005,SPANISH	
L03	3206926-061	21501-001857	LABEL-ADULT VF DOSE, AHA 2005,ITALIAN	
L03	3206926-181	21501-001862	LABEL-ADULT VF DOSE, AHA 2005,SWEDISH	
L03	3206926-201	21501-001863	LABEL-ADULT VF DOSE, AHA 2005,DANISH	
L03	3206926-101	21501-001859	LABEL-ADULT VF DOSE, AHA 2005, DUTCH	
L03	3206926-221	21501-001864	LABEL-ADULT VF DOSE, AHA 2005,FINNISH	
L03	3206926-241	21501-001865	LABEL-ADULT VF DOSE, AHA 2005,NORWEGIAN	
L03	3206926-261	21501-001866	LABEL-ADULT VF DOSE, AHA 2005,POLISH	
L03	3206926-141	21501-001861	LABEL-ADULT VF DOSE, AHA 2005, POR	
L03	3206926-281	21501-001867	LABEL-ADULT VF DOSE, AHA 2005,HUNGARIAN	
L03	3206926-301	21501-001868	LABEL-ADULT VF DOSE, AHA 2005,CZECH	
L03	3206926-321	21501-001869	LABEL-ADULT VF DOSE, AHA 2005,RUSSIAN	
L03	3206926-341	21501-001870	LABEL-ADULT VF DOSE, AHA 2005,CHINESE	
L03	3206926-380	21501-001871	LABEL-ADULT VF DOSE, AHA 05, KOREAN	
L03	3206926-401	21501-001872	LABEL-ADULT VF DOSE, AHA 2005,JAPANESE	

<b>Title</b>	LP20E Service BOM		<b>Page</b>	13 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

<b>A15 Elastomer Keypad - Languages</b>				
<b>Reference Designator</b>	<b>Part Number</b>	<b>Catalog Number</b>	<b>Description</b>	<b>Notes</b>
A07	3200642-106	21300-004598	KEYPAD-ELASTOMER, NO PACING, ENGLISH	
A07	3200642-091	21300-004231	KEYPAD-ELASTOMER, ALL OPTIONS, ENGLISH	
A07	3200642-107	21300-004740	KEYPAD-ELASTOMER, NO PACING, FRENCH	
A07	3200642-092	21300-004755	KEYPAD-ELASTOMER, ALL OPTIONS, FRENCH	
A07	3200642-108	21300-004741	KEYPAD-ELASTOMER, NO PACING, GERMAN	
A07	3200642-093	21300-004712	KEYPAD-ELASTOMER, ALL OPTIONS, GERMAN	
A07	3200642-109	21300-004744	KEYPAD-ELASTOMER, NO PACING, SPANISH	
A07	3200642-094	21300-004713	KEYPAD-ELASTOMER, ALL OPTIONS, SPANISH	
A07	3200642-110	21300-004743	KEYPAD-ELASTOMER, NO PACING, ITALIAN	
A07	3200642-095	21300-004714	KEYPAD-ELASTOMER, ALL OPTIONS, ITALIAN	
A07	3200642-111	21300-004742	KEYPAD-ELASTOMER, NO PACING, SWEDISH	
A07	3200642-096	21300-004715	KEYPAD-ELASTOMER, ALL OPTIONS, SWEDISH	
A07	3200642-112	21300-004748	KEYPAD-ELASTOMER, NO PACING, DANISH	
A07	3200642-097	21300-004716	KEYPAD-ELASTOMER, ALL OPTIONS, DANISH	
A07	3200642-113	21300-004747	KEYPAD-ELASTOMER, NO PACING, DUTCH	
A07	3200642-098	21300-004717	KEYPAD-ELASTOMER, ALL OPTIONS, DUTCH	
A07	3200642-114	21300-004746	KEYPAD-ELASTOMER, NO PACING, FINNISH	
A07	3200642-099	21300-004718	KEYPAD-ELASTOMER, ALL OPTIONS, FINNISH	
A07	3200642-115	21300-004749	KEYPAD-ELASTOMER, NO PACING, NORWEGIAN	
A07	3200642-100	21300-004719	KEYPAD-ELASTOMER, ALL OPTIONS, NORWEGIAN	
A07	3200642-116	21300-004750	KEYPAD-ELASTOMER, NO PACING, POLISH	
A07	3200642-101	21300-004729	KEYPAD-ELASTOMER, ALL OPTIONS, POLISH	
A07	3200642-117	21300-004751	KEYPAD-ELASTOMER, NO PACING, PORTUGUESE	
A07	3200642-102	21300-004720	KEYPAD-ELASTOMER, ALL OPTIONS, PORTUGUESE	
A07	3200642-118	21300-004752	KEYPAD-ELASTOMER, NO PACING, BRAZILIAN	
A07	3200642-103	21300-004721	KEYPAD-ELASTOMER, ALL OPTIONS, BRAZILIAN	
A07	3200642-141	21300-004753	KEYPAD - ELASTOMER, NO PACING, JAPANESE	
A07	3200642-140	21300-008009	KEYPAD - ELASTOMER, ALL OPTIONS, JAPANESE	
A07	3200642-123	21300-004754	KEYPAD-ELASTOMER, NO PACING, CHINESE	
A07	3200642-122	21300-004723	KEYPAD-ELASTOMER, ALL OPTIONS, CHINESE	
A07	3200642-126	21300-006164	KEYPAD-ELASTOMER, LP20 NO PACING, HUN	
A07	3200642-125	21300-006163	KEYPAD-ELASTOMER, LP20 ALL OPTIONS, HUN	
A07	3200642-129	21300-006167	KEYPAD-ELASTOMER, LP20, NO PACING,CZECH	
A07	3200642-128	21300-006166	KEYPAD-ELASTOMER, LP20, ALL OPTIONS, CZECH	
A07	3200642-132	21300-006170	KEYPAD-ELASTOMER, LP20, NO PACING, RUSSIAN	
A07	3200642-131	21300-006169	KEYPAD-ELASTOMER, LP20, ALL OPTIONS,	

<b>Title</b>	LP20E Service BOM		<b>Page</b>	14 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

A15 Elastomer Keypad - Languages				
			RUSSIAN	
A07	3200642-135	21300-006173	KEYPAD-ELASTOMER, LP20, NO PACING, KOREAN	
A07	3200642-134	21300-006172	KEYPAD-ELASTOMER, LP20, ALL OPTIONS, KOREAN	

Standard Paddles Parts List				
Reference Designator	Part Number	Catalog Number	Description	Notes
M35	3200941-004	21330-001024	ADAPTER ASSY-ELECTRODE,HARD PADDLE,PAD PRINTED	
M36	Multiple	Multiple	Charge Button	Refer to Charge Button Languages
M37	3200943-002	21300-008049	BUTTON-SHOCK,HARD PADDLES	
L12	Multiple	Multiple	Paddle Label	Refer to Standard Paddles Label Languages

Standard Paddles Label Languages				
Reference Designator	Part Number	Catalog Number	Description	Notes
L12	3202533-000	21501-000823	LABEL SET-HARD PADDLES,ENGLISH/INTERNATIONAL	
L12	3202533-001	21501-000824	LABEL SET-HARD PADDLES,SPA	
L12	3202533-002	21501-000827	LABEL SET-HARD PADDLES,ITA	
L12	3202533-003	21501-000828	LABEL SET-HARD PADDLES,FIN	
L12	3202533-004	21501-000829	LABEL SET-HARD PADDLES,POL	
L12	3202533-005	21501-000809	LABEL SET-HARD PADDLES,POR/BRA	
L12	3202533-006	21501-000807	LABEL SET-HARD PADDLES,CHI	
L12	3202533-280	21501-001428	LABEL SET-HARD PADDLES, HUNGARIAN	
L12	3202533-300	21501-001429	LABEL SET-HARD PADDLES, CZECH	
L12	3202533-320	21501-001430	LABEL SET-HARD PADDLES, RUSSIAN	
L12	3202533-380	21501-001431	LABEL SET-HARD PADDLES, KOREAN	

<b>Title</b>	LP20E Service BOM		<b>Page</b>	15 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

Charge Button Languages				
Reference Designator	Part Number	Catalog Number	Description	Notes
M36	3200648-015	21300-004886	BUTTON-CHARGE,HARD PADDLES,ENGLISH	
M36	3200648-016	21300-004756	BUTTON-CHARGE,HARD PADDLES,FRENCH	
M36	3200648-017	21300-004757	BUTTON-CHARGE,HARD PADDLES,SPANISH	
M36	3200648-018	21300-004758	BUTTON-CHARGE,HARD PADDLES,ITALIAN	
M36	3200648-019	21300-004759	BUTTON-CHARGE,HARD PADDLES,SWEDISH	
M36	3200648-020	21300-004760	BUTTON-CHARGE,HARD PADDLES,DANISH	
M36	3200648-021	21300-004761	BUTTON-CHARGE,HARD PADDLES,FINNISH	
M36	3200648-022	21300-004762	BUTTON-CHARGE,HARD PADDLES,NORWEGIAN	
M36	3200648-023	21300-004763	BUTTON-CHARGE,HARD PADDLES,POLISH	
M36	3200648-902	21300-004766	BUTTON-CHARGE,HARD PADDLES,GERMAN/DUTCH	
M36	3200648-901	21300-004767	BUTTON-CHARGE,HARD PADDLES,PORTUGUESE,BRAZILIAN	
M36	3200648-024	21300-004764	BUTTON-CHARGE,HARD PADDLES,JAPANESE	
M36	3200648-025	21300-004765	BUTTON-CHARGE,HARD PADDLES,CHINESE	
M36	3200648-280	21300-006328	BUTTON-CHARGE, HARD PADDLES, HUNGARIAN	
M36	3200648-300	21300-006329	BUTTON-CHARGE, HARD PADDLES, CZECH	
M36	3200648-320	21300-006330	BUTTON-CHARGE, HARD PADDLES, RUSSIAN	
M36	3200648-380	21300-006331	BUTTON-CHARGE, HARD PADDLES, KOREAN	

CMM Replacement Parts List				
Reference Designator	Part Number	Catalog Number	Description	Notes
E11	3312422-011	11150-000016	ASSEMBLY - CODE MANAGEMENT MODULE, WI-FI 2.0	
E12	3312422-012	11150-000017	ASSEMBLY-CODE MANAGEMENT MOD,WI-FI AND CO2,2.0	
M38	3313521-001	21300-008092	CODE MANAGEMENT MODULE - BATTERY DOOR	
E13	3311066-002	11141-000162	BATTERY PACK, LITHIUMION, CMM-CODE, MGT, MDL	
E14	3312297-000	21300-008093	MODULE,CAPNOGRAPHY,ORIDION,MICROMEDICO2	
F13	201407-069	21300-000777	SCREW,SELF-SEAL,SELF-LOCK,6-32X0.375	
F14	202253-760	21300-006251	SCREW-M,PH,NYLOK,CS,4-40,.250L	
F15	202253-538	21300-006159	SCREW-LKG,PAN,PHH,4-40,.375,SS,PSVT	
W18	3313521-006	21300-008096	CODE MANAGEMENT MODULE - RF CABLE	
L18	3313521-005	21300-008095	CODE MANAGEMENT MODULE - CO2 MOUNTING SPACER	

<b>Title</b>	LP20E Service BOM		<b>Page</b>	16 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

CMM Replacement Parts List				
M39	3313521-004	21300-008094	CODE MANAGEMENT MODULE - RUBBER GROMMET	
M40	3012120-004	21300-007445	COVER - C02 CONNECTOR	
M28	802885-00	21300-002137	FOOT-MOUNTING,LP8/9/10	
F16	3313521-000	21300-008091	CODE MANAGEMENT MODULE - MOUNTING HARDWARE	

## 5. SERVICE REPAIR KITS

LIFEPAK 20E Service Repair Kits include components that support a particular replacement activity. These repair kits are also listed in 3202718, LP20 Service Replacement Kits. Reference the electronic PLM system Bill of Materials for part numbers and quantity for each assembly.

LP20E Service Repair Kits			
Reference Designator	Part Number	Catalog Number	Description
K01	3202718-000	40402-000001	KIT-REPAIR, TOP CASE HANDLE, LP 20
K02	3202718-001	40402-000002	KIT-REPAIR, FRONT CASE, LP20
K03	3202718-008	40402-000009	KIT REPAIR,DISPLAY FOR UI 3201966-005 AND -006
K04	3202718-009	40402-000010	KIT REPAIR,BACKLIGHT FOR UI 3201966-005 AND -006
K05	3202718-013	40402-000014	KIT-REPAIR, MASIMO SPO2 CONNECTOR, LP20
K06	3202718-021	40402-000017	KIT-REPAIR,PCB STACK,LP20
K07	3202718-024	40402-000020	KIT-REPAIR, CAPACITOR BRACKET, LP20, LP20E
K08	3202718-025	40402-000021	KIT-REPAIR, CAPACITOR, LP20, LP20E
K09	3202718-032	40402-000027	KIT REPAIR,LED DISPLAY FOR UI3201966-007
K10	3202718-041	40402-000034	KIT REPAIR, SPO2 MODULE , LP20, LP20E, ROHS
K11	3202718-042	40402-000035	KIT REPAIR, MNC FLASH SPO2, LP20E, ROHS
K12	3202718-043	40402-000036	KIT REPAIR,THERAPY CONNECTOR, LP20, LP20E,ROHS
K13	3202718-049	40402-000042	KIT REPAIR,THERAPY PCBA,PACING,LP20,4TH ED
K14	3202718-050	40402-000043	KIT REPAIR,THERAPY PCBA, NON-PACING,LP20,4TH ED
K15	3202718-051	40402-000044	KIT REPAIR,OEM PCBA,LP20,4TH ED
K16	3202718-052	40402-000045	KIT REPAIR,PP PCBA,LP20,4TH ED
K17	3202718-057	40402-000050	KIT REPAIR,POST REDUX SC,4TH ED
K18	3202718-058	40402-000051	KIT REPAIR,POST REDUX,PRINTER,4TH ED
K19	3202718-059	40402-000052	KIT REPAIR,POWER MODULE,LP20E,4TH ED V2
K20	3202718-060	40402-000053	KIT REPAIR,POST REDUX UI NO C181,SC,PRINTER,4TH ED
K21	3202718-061	40402-000054	KIT REPAIR,POST REDUX UI NO C181,PRINTER,4TH



<b>Title</b>	LP20E Service BOM		<b>Page</b>	17 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

LP20E Service Repair Kits			
			ED
K22	3202718-062	40402-000055	KIT REPAIR,POST REDUX UI NO C181,4TH ED
K24	3202718-035	40402-000028	KIT REPAIR, THERAPY PCBA, PACING, LP20, ROHS
K25	3202718-036	40402-000029	KIT REPAIR, THERAPY PCBA, NON-PACING,LP20,ROHS
K26	3202718-037	40402-000030	KIT REPAIR, OEM PCBA, LP20, ROHS
K27	3202718-038	40402-000031	KIT REPAIR, PP PCBA, LP20, ROHS
K28	3202718-039	40402-000032	KIT REPAIR, POWER MODULE, LP20E, ROHS
K29	3202718-040	40402-000033	KIT REPAIR, REDUX SC PCBA AND UI PCBA., ROHS
K30	3202718-044	40402-000037	KIT - REPAIR, POST REDUXUI, SC, PRINTER, ROHS
K31	3202718-045	40402-000038	KIT - REPAIR, POST REDUXUI, PRINTER, ROHS
K32	3202718-046	40402-000039	KIT - REPAIR, POST REDUXUI, ROHS
K33	3202718-047	40402-000040	KIT - REPAIR, POST REDUXSC, ROHS
K34	3202718-048	40402-000041	KIT - REPAIR, POST REDUX, PRINTER, ROHS
K35	3202718-054	40402-000047	KIT REPAIR,POST REDUX UI,SC,PRINTER,4TH ED
K36	3202718-055	40402-000048	KIT REPAIR,POST REDUX UI,PRINTER,4TH ED
K37	3202718-056	40402-000049	KIT REPAIR,POST REDUX UI,4TH ED

AED Door/Latch Label Kits				
Reference Designator	Part Number	Catalog Number	Description	Notes
K23	3202360-027	21330-001007	DOOR-AED,W/ LABEL,LP20,ENG/SPA/POR/BRA	AED door (REF 21300-004252) and Latch label (REF 21501-000767)
K23	3202360-028	21330-001008	DOOR-AED,WITH LABEL,LP20,GER/SWE/NOR	AED door (REF 21300-004252) and Latch label (REF 21501-000925)
K23	3202360-029	21330-001009	DOOR-AED,WITH LABEL,LP20,ITA	AED door (REF 21300-004252) and Latch label (REF 21501-000932)
K23	3202360-030	21330-001010	DOOR-AED,WITH LABEL,LP20,FRE/DAN	AED door (REF 21300-004252) and Latch label (REF 21501-000924)
K23	3202360-031	21330-001011	DOOR-AED,WITH LABEL,LP20,DUT	AED door (REF 21300-004252) and Latch label (REF 21501-000930)
K23	3202360-032	21330-001012	DOOR-AED,WITH LABEL,LP20,FIN	AED door (REF 21300-004252) and Latch label (REF 21501-000929)
K23	3202360-033	21330-001013	DOOR-AED,WITH LABEL,LP20,POL	AED door (REF 21300-004252) and Latch label (REF

<b>Title</b>	LP20E Service BOM		<b>Page</b>	18 of 18
<b>Doc Type</b>	Service Product Drawing	<b>Doc ID</b>	3343496	<b>Doc Rev</b> AC
<b>Doc State</b>	Released	<b>State Date</b>	9/1/2023	

AED Door/Latch Label Kits				
Reference Designator	Part Number	Catalog Number	Description	Notes
				21501-000928)
K23	3202360-034	21330-001014	DOOR-AED,WITH LABEL,LP20,CHI	AED door (REF 21300-004252) and Latch label (REF 21501-000927)
K23	3202360-035	21330-001015	DOOR-AED,WITH LABEL,LP20,JAP	AED door (REF 21300-004252) and Latch label (REF 21501-000926)