



Mattress Assembly - 2235-100-001



Part		Description	Dimensions
A	2235-100-001	Main Mattress	36" x 81" x 6.5"
B	2235-100-004	Foot Extension	8" x 35.5" x 6.5"
C	2235-100-002	Right Side Extension	6" x 81" x 6.5"
D	2235-100-003	Left Side Extension	6" x 81" x 6.5"
E	2235-100-009	Right Corner Extension	6" x 8" x 6.5"
F	2235-100-010	Left Corner Extension	6" x 8" x 6.5"